lumie[®] vitamin L[™]

Safety

WARNING - To reduce the risk of fire, electric shock or injury:

- If the unit is damaged contact Lumie before use.
- Keep away from water and damp. 
- For indoor use only. 
- The unit gets warm when the light is on, so make sure that it is not covered.
- Use only the mains power adaptor supplied.
- Keep the mains power adaptor away from heated surfaces. If it is damaged it should be replaced with an approved Lumie product.
- To power down or reset the unit, remove the mains power adaptor from the wall socket.
- Do not attempt to service or modify any part of the unit.
- Do not shorten or cut the power cord.
- If the unit has been stored outside the recommended temperature range (5 °C to 35 °C / 41 °F to 95 °F), allow it to stand for at least 1 hour before applying power to avoid condensation or overheating.
- Not to be used by or on individuals who are unconscious, anaesthetized or incapacitated e.g. under the influence of alcohol.

CAUTION

- Do not drop the unit on the floor, hit it hard or expose it to other shocks.
- Place the unit on a stable, level and non-slippery surface.
- Always unplug the unit before cleaning and allow it to cool down.
- Not to be used with light dimmers, timers, motion detectors, voltage transformers or extension cords.
- Intended for household use, including similar use in hotels. Do not use attachments not recommended by Lumie.

Electromagnetic compatibility

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

For further detailed information on the electromagnetic environment in which the device can be operated refer to the website: www.lumie.com/pages/e-m-c.

Responsible Party Address:

Violette Engineering Corporation
6731 Whittier Avenue
McLean, VA 22101
info@violettecorp.com

Disposal

At the end of its life, do not throw the unit out with normal domestic waste. Please take it to an official recycling point to be disposed of responsibly.

	Serial number		For indoor use only		Standby
	Manufacturer: Outside In (Cambridge) Ltd, 3 The Links, Bar Hill, Cambridge, CB23 8UD, UK.		Read instructions for use		Recyclable
	Keep dry		Double-insulated (Class II electrical)		

Precautions

Research shows that light therapy is very safe. However, please check with your doctor first if any of the following apply:

- You have a major eye problem or have had eye surgery;
- You have been advised to avoid bright light or are taking medication that may cause photosensitivity;
- You have been diagnosed with depression or other mood/sleep disorders. This is especially important if you are also taking medication for this, since light therapy can also affect your mood.

A few people experience side effects but these are generally mild (including headaches, eye strain, nausea and hyperactivity) and can usually be solved simply by switching off the light. Over the next few days, gradually build up the time you spend in front of it to find a treatment time that works best for you.

Some people complain of difficulty sleeping and often this is down to using the light too close to bedtime.

If you have any problems while using Vitamin L, or if you've been using the light for a week but haven't noticed any benefit, please contact Lumie or speak to your doctor.

IMPORTANT: the guidelines in this leaflet are not medical directions for treating any condition.

Introduction

Using Vitamin L will provide you with the bright light your body needs during the dark winter months. It can help to boost your energy levels, put you in a better mood and make you feel more awake¹. Light therapy is safe and is proven to help alleviate the winter blues.


Please read these instructions carefully to get the most out of your Vitamin L. 

1. Acta Psychiatr Scand. 2001 Apr;103(4):267-74.

Disclaimer: 'Vitamin L' is simply a product name where 'L' refers to 'light'. The product has no connection to any real vitamin and Lumie does not claim any effects on vitamin levels in the body as a result of using this product.

Using Vitamin L

- Connect the mains power adaptor to the socket at the back of the unit.
- Plug in the unit.
- Press  to turn on/off.
- Slot the stand into the back of the unit. For extra stability, reverse the stand when switching between portrait and landscape use:

- Position Vitamin L at arm's length (about 50 cm/20 in) so the light is shining towards your face.

Although you don't need to stare at it, the light must reach your eyes to have any effect. This means you can use Vitamin L whilst you're working at the computer, watching TV, reading, etc. Dark glasses, tinted lenses or closing your eyes will all limit the effects of light therapy.

You should notice positive benefits after using Vitamin L for three or four days in a row.

When should I use it?

It's best to tackle the winter blues as soon as possible. Look out for early signs - often in September or October - and if you start to feel lethargic, anxious or irritable, start using Vitamin L.

Most people can skip a day here and there, especially when the weather is bright. As spring comes around you will probably start to cut back your light therapy; your body will tell you if it's too soon and if you notice your symptoms returning simply use your Vitamin L again for a few days.

Light therapy is proven to put you in a better mood, to boost alertness and to help you feel more awake. Use Vitamin L wherever and whenever you need to restore concentration and give yourself a mental boost.

What time of day?

If you're using Vitamin L to help with the winter blues, bright light in the morning has been shown to be most effective so try that first. For some people, early evening light exposure is beneficial.

Whether you're using Vitamin L for winter blues or as a general pick-me-up, once you've found a pattern that's right for you, try to use it around the same time every day.

I find it difficult to get going in the morning

Switch on as soon as you can. If you don't have much time in the morning, then you should aim for a top-up or full session in the afternoon.

I feel sleepy or fall asleep too early

Use Vitamin L in the late afternoon/early evening. It's best to avoid bright light too close to bedtime - within about three hours - as you may find it difficult to get to sleep.

How long?

Everyone's different, so position Vitamin L at arm's length (about 50 cm/20 in) and start with 30 minutes a day. Some people will find this is enough for effective treatment. If you need to, increase your exposure each day until you find a treatment time that works for you.

Distance (cm/in)	Light received (lux)	Treatment time (minutes)
20/6.3	10,000	30
30/11.8	5,000	30–60
50/19.7	2,500	60–120

You can get your daily dose of Vitamin L all in one go or use it on-and-off throughout the day for a cumulative effect. If you're more comfortable with Vitamin L further away, then less light will reach your eyes and you'll need to use it for longer.

There's no harm in using the light for longer than suggested. If you find your eyes ache slightly, try using it further away or switch it off for a while.

LED lifetime

The LEDs in this product have a typical life of 30,000 hours. Although it's unlikely, it's possible for individual LEDs to fail. This will make only a negligible difference to light output and you can continue to use the product.

Maintenance

Clean only with a dry or lightly moistened cloth.

Replacement mains power adaptors are available from Lumie.

Warranty

This product is covered under warranty against any manufacturing defect for 3 years from the date of purchase. This is in addition to your statutory rights. If you need to return the product for servicing, use the original packaging (including internal packing) if possible. Make sure the mains power adaptor is prevented from moving in transit and damaging the unit. If the unit is received damaged, through misuse or accidental damage, then we may have to charge to repair it; we will contact you first. Please contact Lumie for further advice.

Technical description

Subject to change without notice.

Lumie Vitamin L

100-240 Vac, 50/60 Hz, 0.5 A max power adaptor (supplied) 

Operating temperature range: 5 °C to 35 °C / 41 °F to 95 °F

Storage/transport temperature range: -25 °C to +70 °C / -13 °F to +158 °F

Contact us

If you have any suggestions that you feel could help us improve our products or service, or if there are problems you think we should know about, please call the Lumie Careline on +44 (0)1954 780500 or email info@lumie.com.

Lumie, 3 The Links, Trafalgar Way, Bar Hill, Cambridge, CB23 8UD, UK.

Lumie is a registered trademark of:

Outside In (Cambridge) Limited, incorporated in England and Wales.

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© Lumie 2019
Outside In (Cambridge) Limited, constituida en Inglaterra y Gales.