lumie®

brazil™

BUL (Bulgarian)

CZE (Czech)

DAN (Danish)

DEU (German)

DUT (Dutch)

ENG (English)

EST (Estonian)

FIN (Finnish)

FRA (French)

GRE (Greek)

HRV (Croatian)

HUN (Hungarian)

ITA (Italian)

LAV (Latvian)

LIT (Lithuanian)

NOR (Norwegian)

POL (Polish)

POR (Portuguese)

RUM (Romanian)


SLO (Slovak)

SLV (Slovenian)

SPA (Spanish)

SWE (Swedish)

TUR (Turkish)

lumie® brazil™

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
GB 880 9837 71 (UK)

NL 826301381B01 (EU)

 Advena Ltd, Tower Business Centre, 2nd Floor, Tower Street, Swatar, BKR 4013 Malta.

lumie® brazil™

Safety

If the unit is damaged, contact Lumie before use.

Keep away from water and damp. 

For indoor use only. 

The unit gets warm when the light is on, so make sure that it is not covered.

Keep the cable away from heated surfaces. If the cable is damaged it must be replaced by Lumie, a Lumie distributor or suitably qualified person.

The only way to disconnect the unit from the power supply is to remove the plug from the mains socket. You should therefore ensure the socket is easily accessible.

Only to be used by adults and children over 7 years old. Children using the light, or playing near it, should be supervised.

WARNING: If the wall socket used to power the device has poor connections, the plug becomes hot. Make sure you plug the device into a properly installed wall socket to avoid the risk of fire and burns.

Do not attempt to service or modify any part of the unit.

If the unit has been stored outside the recommended operating temperature range (5 °C to 35 °C), allow it to stand for at least 1 hour before applying power to avoid condensation or overheating.

Do not operate the device in the presence of flammable gases or if there is a high concentration of oxygen.

CAUTION: Do not subject the device to heavy shocks to avoid risk of damage to the lamp.

Not to be used by or on individuals who are unconscious, anaesthetised or incapacitated e.g. under the influence of alcohol.



Electromagnetic compatibility

This product conforms to the electromagnetic compatibility standard for medical devices and is unlikely to cause interference affecting other electrical and electronic products in its vicinity. However, radio transmitting equipment, mobile phones and other electrical/electronic equipment can interfere with this product; if this happens simply move the product away from the source of interference.

For further detailed information on the electromagnetic environment in which the device can be operated refer to the website: www.lumie.com/pages/e-m-c.

Disposal

Do not throw out the bulbs or the unit with normal domestic waste. Take them to an official recycling point to be disposed of responsibly. 

	Serial number		Read instructions for use		Double-insulated (Class II electrical)
	European Declaration of Conformity: Medical Device Class IIa		Keep dry		AC (alternating current) power supply
	Medical Device		For indoor use only		Push ON/OFF
	Unique Device Identifier		Caution		Fuse
	Manufacturer		Do not place in household waste. Dispose of at a recycling centre.		Storage/transport relative humidity
	Authorised representative in the European Community		Member of the 'Green Dot' recycling scheme		Storage/transport temperature
	Importer		Recyclable		UK Conformity Assessed

Precautions

Research shows that light therapy is very safe. However, please check with your doctor first if any of the following apply:

- You have a major eye problem or have had eye surgery;
- You have been advised to avoid bright light or are taking medication that may cause photosensitivity;
- You have been diagnosed with depression or other mood/sleep disorders. This is especially important if you are also taking medication for this, since light therapy can also affect your mood.

A few people experience side effects but these are generally mild (including headaches, eye strain, nausea and hyperactivity) and can usually be solved simply by switching off the light. Over the next few days, gradually build up the time you spend in front of it to find a treatment time that works best for you.

Some people complain of difficulty sleeping and often this is down to using the light too close to bedtime.

If you have any problems while using Brazil, or if you've been using the light for a week but haven't noticed any benefit, please contact Lumie or speak to your doctor.

IMPORTANT: the guidelines in this leaflet are not medical directions for treating any condition. If your doctor has suggested using light therapy, follow their instructions and discuss any effects with them.

Introduction

Using Brazil will provide you with the bright light your body needs during the dark winter months. It can help to boost your energy levels, put you in a better mood and make you feel more awake¹. Light therapy is safe and is proven to help alleviate the winter blues and SAD (Seasonal Affective Disorder)².

Please read these instructions carefully to get the most out of your Brazil. 

1. Acta Psychiatr Scand. 2001 Apr;103(4):267-74.
2. Am J Psychiatry. 2005 Apr;162(4):656-62.

Using Brazil

- Remove the packing around the bulbs before first use. See *Troubleshooting*, page 7, for how to remove the screen.
- Plug in the unit and switch on at the mains.
- Press  to turn on/off.
- Position Brazil at arm's length (about 50 cm) so the light is shining towards your face.

Although you don't need to stare at it, the light must reach your eyes to have any effect. This means you can use Brazil whilst you're working at the computer, watching TV, reading, etc. Dark glasses, tinted lenses or closing your eyes will all limit the effects of light therapy.

You should notice positive benefits after using Brazil for three or four days in a row.

When should I use it?

It's best to tackle the winter blues and SAD as soon as possible. Look out for early signs - often in September or October - and if you start to feel lethargic, anxious or irritable, start using Brazil.

Most people can skip a day here and there, especially when the weather is bright. As spring comes around you will probably start to cut back your light therapy; your body will tell you if it's too soon and if you notice your symptoms returning simply use your Brazil again for a few days.

Light therapy is proven to put you in a better mood, to boost alertness and to help you feel more awake. Use Brazil wherever and whenever you need to restore concentration and give yourself a mental boost.

What time of day?

If you're using Brazil to help with the winter blues or SAD, bright light in the morning has been shown to be most effective so try that first. For some people, early evening light exposure is beneficial.

Whether you're using Brazil for winter blues, SAD or as a general pick-me-up, once you've found a pattern that's right for you, try to use it around the same time every day.

I find it difficult to get going in the morning

Switch on Brazil as soon as you can. If you don't have much time in the morning, then you should aim for a top-up or full session in the afternoon.

I feel sleepy or fall asleep too early

Use Brazil in the late afternoon/early evening. It's best to avoid bright light too close to bedtime - within about three hours - as you may find it difficult to get to sleep.

How long?

Everyone's different, so position Brazil at arm's length (about 50 cm) and start with 30 minutes a day. Some people will find this is enough for effective SAD treatment. If you need to, increase your exposure each day until you find a treatment time that works for you.

Distance (cm)	Light received (lux)	Treatment time (minutes)
35	10 000	30
50	5 000	30–60
80	2 500	60–120

You can get your daily dose of light therapy all in one go or use it on-and-off throughout the day for a cumulative effect. If you're more comfortable with Brazil further away, then less light will reach your eyes and you'll need to use it for longer.

There's no harm in using the light for longer than suggested. If you find your eyes ache slightly, try using it further away or switch it off for a while.

Maintenance

Clean only with a dry or lightly moistened cloth. Always unplug the unit first and let it cool if necessary.

The bulbs will slowly lose their brightness so, to maintain maximum output and benefit, we recommend you replace them every two to three years. Replacement bulbs are available from Lumie.

Troubleshooting

The light is less bright or doesn't come on at all

- One or more of the bulbs may have failed. Switch off and unplug the unit. Let it cool if necessary.
- Insert a coin or something similar into the slot at the top centre of the plastic screen and prise it away from the case. Detach from the sides and unhook the screen at the bottom.
- Check that the bulbs are firmly in the sockets. If you need to adjust them, plug in the unit and switch on afterwards to see if this has worked.
- To replace the screen, hook the bottom edge first and push the top and sides to click it back into position.

The bulbs are firmly in place but it hasn't fixed the problem

- If none of the bulbs are working, check the fuses (see below). If the fuses are OK, then you may need to replace the bulbs.
- To remove a bulb, push the red button down hard and slide the bulb out of the socket.
- Gently slide the replacement bulb into the socket until it clicks into position.

I've replaced the bulbs but it's still not working

- Check the internal fuses on each side of the on/off switch - it will be obvious if either has blown. Replace with a 5x20mm cartridge fuse, type F3.15AL250V.
- It is also possible that the plug fuse might have blown. Replace with a 3A fuse to BS1362.

If this doesn't work or if you have any other problems with your Brazil, please call Lumie.

Note: Any serious incident relating to the unit should be reported to Lumie and to the competent authority in the country where you live.

Warranty

This unit is covered under warranty against any manufacturing defect for 3 years from the date of purchase. This is in addition to your statutory rights. If you need to return the product for servicing, use the original packaging (including internal packing) if possible. Make sure the mains power adaptor is prevented from moving in transit and damaging the unit. If the unit is received damaged, through misuse or accidental damage, then we may have to charge to repair it; we will contact you first. Please contact Lumie for further advice.

Technical description

Subject to change without notice.

Lumie Brazil

230 Vac, 50 Hz, 125 W 

Lamps: 3 x PLL/36 W/840

Operating conditions:

Temperature: 5 °C to 35 °C

Relative humidity: 15% to 93% (non-condensing)

Atmospheric pressure: 70 kPa to 106 kPa

Storage/transport conditions:

Temperature: -25 °C to +70 °C

Relative humidity: 15% to 93% (non-condensing)

Atmospheric pressure: 25 kPa to 106 kPa

CE

European Declaration of Conformity: Medical Device Class IIa

**UK
CA**

0120 UK Conformity Assessed


MD

Medical Device

Optical Radiation Risk Group: Group 1

Maximum Optical Radiation Output: 10 000 lux + 30% at 35 cm

Relative Spectral Irradiance as shown below:

Contact us

If you have any questions or problems regarding this product, or have suggestions that you feel could help us improve our products or service please email info@lumie.com.

 Outside In (Cambridge) Limited, 3 The Links, Trafalgar Way, Bar Hill, Cambridge, CB23 8UD, UK.

Lumie is a registered trademark of:


Outside In (Cambridge) Limited, incorporated in England and Wales.

Registered Number: 2647359.

VAT No:

GB 880 9837 71 (UK)

NL 826301381B01 (EU)

 Advena Ltd, Tower Business Centre, 2nd Floor, Tower Street, Swatar, BKR 4013 Malta.

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© Lumie 2025
Outside In (Cambridge) Limited, 3 The Links, Trafalgar Way, Bar Hill, Cambridge, CB23 8UD, Birleşik Krallık.
İngiltere ve Galler'de kurulmuş olan Outside In (Cambridge) Limited.
NL 826301381B01 (AB)
Outside In (Cambridge) Ltd