lumie[®] bodyclock[™] rise 100

Safety

If the unit is damaged contact Lumie before use.

Keep away from water and damp.

For indoor use only.

The unit gets warm when the light is on, so make sure that it is not covered.

Use only the mains power adaptor supplied.

Keep the mains power adaptor away from heated surfaces. If either are damaged they should be replaced with an approved Lumie product.

To power down or reset the unit, remove the mains power adaptor from the wall socket.

Do not attempt to service or modify any part of the unit.

If the unit has been stored outside the recommended temperature range (5 °C to 35 °C), allow it to stand for at least 1 hour before applying power to avoid condensation or overheating.

Not to be used by or on individuals who are unconscious, anaesthetised or incapacitated e.g. under the influence of alcohol.

Electromagnetic compatibility

This product conforms to the electromagnetic compatibility (EMC) standard for medical devices and is unlikely to cause interference affecting other electrical and electronic products in its vicinity. However, radio transmitting equipment, mobile phones and other electrical/electronic equipment can interfere with this product; if this happens simply move the product away from the source of interference. In rare cases where product function is lost due to EMC effects, disconnect the product from power source for 5 seconds and then restart. If the problem persists, please contact Lumie.

For further detailed information on the electromagnetic environment in which the device can be operated refer to the website: www.lumie.com/pages/e-m-c.

Disposal

At the end of its life, do not throw the unit out with normal domestic waste. Please take it to an official recycling point to be disposed of responsibly.

 SN	Serial number		For indoor use only		Double-insulated (Class II electrical)
	European Declaration of Conformity: Medical Device Class I		Do not place in household waste. Dispose of at a recycling centre.		DC (direct current) power supply
	Manufacturer: Outside In (Cambridge) Ltd, 3 The Links, Bar Hill, Cambridge, CB23 8UD, UK.		Read instructions for use		Member of the 'Green Dot' recycling scheme
	Keep dry		Caution		Recyclable

Introduction

Lumie Bodyclock wakes you naturally with a sunrise. Waking with light has been shown to improve mood, energy, productivity and the quality of sleep and awakening.¹

Wake-up lights can also help sufferers of SAD (Seasonal Affective Disorder) and winter blues cope with dark winter mornings.²

At bedtime, a fading sunset helps you naturally unwind by promoting the production of the sleep hormone, melatonin.³

Use Lumie Bodyclock all year round to encourage a regular sleep-wake pattern.

¹ Acta Psychiatr Scand. 1993 Jul;88(1):67-71.

² J Affect Disord. 2015 Jul 15;180:87-9.

³ Sleep Biol Rhythms. 2016 Jan;14(1):47-53.

Getting started

The ideal place for Lumie Bodyclock is on a bedside table near a mains socket.

- Connect the mains power adaptor to the socket in the base of the unit.
- Plug into a mains power supply and switch on.

The display will flash until the current time is set:

Program

The hours will flash.

- Use  or  to set the hours.
- Press  to confirm and continue.

The minutes will flash.

- Use  or  to set the minutes.
 - Press  to confirm and exit to the current time.
-

Bodyclock Rise 100 is preset with a 07:00 alarm time, 30-minute sunrise and 30-minute sunset.

To change the alarm time or any other settings, see *Alarm time and sunrise options*, page 5, and *Sunset options*, page 7.

Alarm time and sunrise options

- Press and hold .

A symbol on the display shows which option you're setting. Continue through the following sequence to set the alarm time and sunrise options:

Alarm time

- Use  or  to set the hours.
- Press .
- Use  or  to set the minutes.
- Press  to confirm and continue.

Sunrise light intensity

- Use  or  to set the final brightness of the sunrise: 0 for no sunrise or a level from 1 to 5.
- Press  to confirm and continue.

Alarm beep

- Use  or  to set the volume of the optional alarm beep: 0 for no alarm beep or a level from 1 to 5.

If you set the sunrise light intensity to 0 you have to set an alarm beep instead.

- Press  to confirm and exit to the current time.
-

Using the sunrise alarm

You must activate the alarm before you go to sleep each time.

To wake with a sunrise alarm

- Press .

The display confirms the alarm time and shows  to confirm the alarm is activated; for example:

You'll wake with a gradual sunrise and any other sunrise options you have set. See also *Alarm time and sunrise options*, page 5.

To deactivate the sunrise alarm

- Press .

 disappears from the display.

To snooze

- Tap the top of the lamp cover.

The sunrise, and alarm beep if set, will repeat after 9 minutes.

To turn off the sunrise and alarm beep

- Press any button.

 disappears from the display.

- Press and hold  to turn off the light.

Timeout: If the light and alarm beep are not turned off manually, they will turn off automatically after 90 minutes.

Sunset options

- Press and hold .

A symbol on the display shows which option you're setting. Continue through the following sequence to set the sunset options:

Nightlight

- Use  or  to set a nightlight: 0 for no nightlight (sunset fades to off) or a level from 1 to 5 (sunset fades to nightlight).
- Press  to confirm and exit to the current time.


Using the fading sunset

When you're ready to sleep, a fading sunset can help you to relax and unwind. See also *Sunset options* above.

To activate the sunset

- Use  or  to set the light level.
- Press .

The display shows  to confirm the sunset is activated:

The light will start to fade and any other sunset options you have set will be activated.

To deactivate the sunset

- Press .


 disappears from the display.

Display brightness

The display automatically adjusts to match the light level in your room and turns OFF when the room is dark and you're ready to sleep.

You can change this setting if you prefer the display to stay on at LOW brightness when the room is dark:

- Use the switch in the base of the unit to set OFF or LOW:

To check the time when the display is off

- Tap the top of the lamp cover.

The display will change to low brightness for a few seconds.

Changing the current time

Bodyclock Rise 100 will keep time for 30 minutes without power. All other settings are permanently saved.

- Press and hold  until  appears.

Continue through the following sequence to set the time:

	Program	The hours will flash. <ul style="list-style-type: none">• Use  or  to set the hours.• Press  to confirm and continue. The minutes will flash. <ul style="list-style-type: none">• Use  or  to set the minutes.• Press  to confirm and exit to the current time.
--	---------	---

Bedside light

Lumie Bodyclock can be used as a fully dimmable bedside light.

- Use  or  to set the light level.
- Press and hold  to turn off the light.

Timeout: If the bedside light is not turned off manually, the light will turn off automatically after 12 hours.

Maintenance

- Clean only with a dry or lightly moistened cloth.
- Replacement mains power adaptors are available from Lumie.

To restore the default factory settings

- Turn off the power at the mains.
- Press and hold  and turn the power back on.
- Release the button when the software version number shows.

LED lifetime

The LEDs in this product have a typical life of 30,000 hours. Although it's unlikely, it's possible for individual LEDs to fail. This will make only a negligible difference to light output and you can continue to use the product.

Warranty

This unit is covered under warranty against any manufacturing defect for 3 years from the date of purchase. This is in addition to your statutory rights. If you need to return the product for servicing, use the original packaging (including internal packing) if possible. Make sure the mains power adaptor is prevented from moving in transit and damaging the unit. If the unit is received damaged, through misuse or accidental damage, then we may have to charge to repair it; we will contact you first. Please contact Lumie for further advice.

Technical description

Subject to change without notice.

Bodyclock Rise 100

Mains power adaptor (supplied): 

Input: 100–240 Vac, 50/60 Hz, 0.4 A max.; Output: 12 V \equiv 1 A (12 W)

Operating temperature range: 5 °C to 35 °C

Storage/transport temperature range: -25 °C to +70 °C

Medical Device Class I 

Contact us

If you have any suggestions that you feel could help us improve our products or service, or if there are problems you think we should know about, please call the Lumie Careline on +44 (0) 1954 780500 or email info@lumie.com.

Lumie, 3 The Links, Trafalgar Way, Bar Hill, Cambridge, CB23 8UD, UK.

Lumie is a registered trademark of:

Outside In (Cambridge) Limited, incorporated in England and Wales.

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© Lumie 2019
Outside In (Cambridge) Limited, constituida en Inglaterra y Gales.